


# Scan Me to Sign-In!





MODERN MOBILITY

## **Unlocking Independence:**

**SMART Technology Solutions for Wheelchair Users with Visual Impairments**




## Authored and Presented By

**Lindsey Sharpe, PT, DPT, ATP**  
LUCI: Regional Director of Clinical  
Education  
lindsey@luci.com

# Learning Objectives


- ✓ Participants will identify at least 2 common safety and accessibility challenges experienced by wheelchair users with combined visual and mobility impairments.
- ✓ Participants will analyze at least 2 research studies that examine both the benefits of self-generated mobility and the risks of limited access to power mobility for individuals with visual and motor impairments.
- ✓ Participants will describe at least 2 features of current SMART wheelchair technologies and their role in promoting safe, independent wheelchair use.
- ✓ Participants will apply at least 2 clinical strategies for evaluating and recommending technology features that enhance confidence, safety, and independence in power wheelchair users with visual impairments.



# Overview for Today



The Need



The Technology



The Person

# The Need

Power Wheelchair Limitations: Efficiency and Safety

# Main Goal of Power Wheelchairs

Designed to provide an alternative means of moving throughout the environment for people who are unable to safely and/or efficiently **AMBULATE** or **SELF-PROPEL** an optimally configured manual wheelchair

**GOAL:**


Access to self-generated mobility!



# Modes of Self-Generated Mobility

At each level, ask: Is this mobility mode functional, efficient, and safe?

➤ Can the person manage over 24 hr period? Can the person navigate school/ALF?



# Research Shows


## The Value of Independence is Immeasurable

<p><b>Greater Independence</b></p>	<p>PMD use is associated with an increased frequency of grocery shopping and going for “walks”, and an increased frequency of instrumental activities of daily living, such as going to a restaurant, posting letters, going to the bank, and visiting family and friends.</p> <p><i>Source: Pellichero, A., et al. Relationships between Cognitive Functioning and Powered Mobility Device Use: A Scoping Review. Int. J. Environ. Res. Public Health 2021.</i></p>
<p><b>Increased Brain Development</b></p>	<p>In children, PMD use contributed to the development of cognitive and play skills while increasing independence and social interactions.</p> <p><i>Source: Pellichero, A., et al. Relationships between Cognitive Functioning and Powered Mobility Device Use: A Scoping Review. Int. J. Environ. Res. Public Health 2021.</i></p>
<p><b>Promotes Self-Reliance</b></p>	<p>Independent mobility increases vocational and educational opportunities, reduces dependence on caregivers and family members, and promotes feelings of self-reliance. Reductions in functional mobility are linked with reduced participation and loss of social connections...decreases in mobility can lead to feelings of emotional loss, reduced self-esteem, isolation, stress, and fear of abandonment</p> <p><i>Source: Simpson R, et al. .How many people would benefit from a smart wheelchair? Journal of Rehabilitation Research and Development. 2008</i></p>
<p><b>It is a Human Right</b></p>	<p>According to a social justice lens, each and every individual, regardless of disability, has a fundamental right to self-directed mobility in order to fully participate in life as defined by the International Classification of Functioning, Disability and Health framework.</p> <p><i>Source: Samuel W et al. Factors predicting attitudes toward self-directed mobility. Disability and Health Journal. 2018</i></p>

# Real Talk: Driving a Power Wheelchair is Challenging

**Wheelchair drivers, their families, caregivers + clinicians** report reluctance to pursue power wheelchairs because:


- They worry that the driver won't be safe
- They worry that the driver will hurt someone else
- They worry that the driver will damage the environment, i.e. walls and doorframes in their home




“Often, when medical professionals choose not to prescribe wheelchairs for their patients, they have not taken into account the severe limitations put on independent travel when a person cannot drive or walk more than a block or two. Invite the PT and OT to attend a few lessons; it may change their view on the need for a powered aide.”

American Printing House for the Blind. Orientation & Mobility for Wheelchair Users: Chapter 01. Retrieved January 13, 2026, from <https://tech.aph.org/omwc/xhtml/chapter-01.xhtml>.

## O+M Assessment:

- 
- Intake Interview
  - The Wheelchair
  - Movement
  - Turns Around Obstacles
  - Emergency Stops
  - Speed
  - Turns at Corners
  - Doors
- Elevators
  - Building Ramps
  - Sidewalks
  - Curbs and Curb Ramps
  - Street Crossings
  - Transportation
  - Route Planning
  - Feasibility of Using Powered Chairs
  - Consultations

Now, we layer on visual impairments and potential need for a white cane



Operating a power wheelchair is a lot like driving a car...but in MUCH more challenging environments!

What can get in the way of SAFE driving?

- **Visual Impairments**
- **Poor Motor Control**
- **Slowed Reaction time**
- **Highly Distractible Environments**

“It is RESNA’ s position that age, **limited vision** or cognition, behavioral issues, and the ability to walk or propel a manual wheelchair short distances should not, in and of themselves, be used as discriminatory factors against providing powered mobility for children.”

*Source: Rosen L, et al. RESNA position on the application of power mobility devices for pediatric users. Assist Technol. 2023*



# Challenges from the Built Environment

## Transportation

- Heavy → requires accessible transportation
- Narrow ramps
- Locking systems in vehicles

## Accessibility

- Getting in and out of buildings and navigating throughout tight and/or crowded spaces → including home



## Challenges with White Cane + Manual Wheeled Mobility

- For a person needing wheeled mobility, manual wheelchair propulsion while using a white cane becomes more complex
  - e.g., one hand on wheels, turning, wheel-chair momentum/stopping, etc.



# Challenges from the Power Wheelchair

## Inefficiency

- Maintaining a 'straight' path can be difficult, especially with alternative drive controls
- Caster flips can divert the chair



*Tracking technologies can help but are not standard. (And not always requested!)*

## Safety Concerns

- The chair design cannot prevent:
  - Collisions
  - Detect and evade a drop-off
  - Tip over due to a steep angle



*Reliance is solely on the driver to note and avoid potential hazards.*

## Consider the Clients You Work With...

10-40% denied access to power mobility...


(Fehr L, et al. 2000)

- ✓ Sensory impairments
- ✓ Poor motor function
- ✓ Cognitive deficits



40% of existing power wheelchair users are struggling

(Fehr L, et al. 2000)



---

# Combined Visual + Mobility Deficits

---

**The topic of teaching Orientation and Mobility (O&M) to people with orthopedic impairments has received little attention in blindness-related literature.**

*American Printing House for the Blind. Orientation & Mobility for Wheelchair Users: Chapter 01. Retrieved January 13, 2026, from <https://tech.aph.org/omwc/xhtml/chapter-01.xhtml>.*

---

**Almost 10% of all individuals who are legally blind also have a mobility impairment.**

*Simpson R, et al. A prototype power assist wheelchair that provides for obstacle detection and avoidance for those with visual impairments. J Neuroeng Rehabil. 2005.*

---

**...visual impairments can make obstacle detection, spatial navigation, hand-eye coordination and use of orientation cues difficult for wheelchair users.**

*Simpson R, et al. .How many people would benefit from a smart wheelchair? Journal of Rehabilitation Research and Development, 2008.*

---

**...many visually impaired wheelchair users are dependent on others for mobility.**


*Simpson R, et al. A prototype power assist wheelchair that provides for obstacle detection and avoidance for those with visual impairments. J Neuroeng Rehabil. 2005.*

---

**Forces of impact from tip and roll accidents result in significant risk for mild to severe head injury, depending on chair position and restraint at the time of incident**

*Source: Brett Erickson, et al. The dynamics of electric powered wheelchair sideways tips and falls: experimental and computational analysis of impact forces and injury. Journal of NeuroEngineering and Rehabilitation, 13(1):20, 2016.*

## Gaps in the Existing Research



# Solution?

## SMART Technology

When combined with SMART technology solutions, wheelchairs can provide increased independence, efficiency and safety by protecting the driver, others around them, and the environment.



# The Technology

Definitions, Terminology, and Classification

# What makes a mobility device SMART?



**SELF-MONITORING:** Sensor that can provide environmental surveillance and detect, for example, obstacles, drop-offs, or inclines.



**ANALYSIS:** Using the data generated by the sensors, a SMART solution analyzes this information according to the user's customized preferences.



**REPORTING:** This analyzed data is then interpreted and “reported” back to the user so they themselves/or in combination with their mobility device can adjust their pathway.



**TECHNOLOGY:** The technology includes both the hardware (the sensors on the chair!) and the software, which provides the analysis and reporting functions.

## SMART in the Context Power Wheelchairs

*A Smart Wheelchair is integrated or retrofitted self-monitoring technology for a power wheelchair that provides enhanced, independent **mobility** to a wheelchair user, can collect and report **user health and wellness data** and provides **connectivity** to integrate with the connected world.*

-Michelle L. Lange, OTR/L, ABDA, ATP/SMS

# SMART WHEELCHAIR (SWC) CONTINUUM

## SWC Level 0

Warning Systems

No Automation/  
Intervention

## SWC Level 1

Driver Assistance

Single function assisted navigation

## SWC Level 2

Advanced Driver Assistance

Multiple function assisted navigation

## SWC Level 3

Conditional Automation

Autonomously navigate through a specific process/under specific conditions

## SWC Level 4

Highly Autonomous System

The wheelchair is fully autonomous for an entire trip in specific environments

## SWC Level 5

Fully Autonomous System

The wheelchair can navigate without human input in all environments

# SMART WHEELCHAIR (SWC) CONTINUUM

SWC Level	Name	Definition	Human Role	Functional Example	Product Feature Example	Product Example
<b>Human driver monitors the driving environment</b>						
0	<b>Warning Systems:</b> no automation/intervention	A warning system that monitors, and alerts, or provides additional feedback to the driver of potential hazards but does not affect user drive inputs. The system does NOT intervene but relies solely on the driver	Operator-the user is always in complete control.	The system does NOT intervene- reliance is solely on the driver to respond appropriately/timely to the warning given	Sensors that warn driver of potential collisions or other hazards	<p><b>Braze Mobility:</b> collision warning systems (auditory, visual, and haptic warnings)</p> <p><b>ASL:</b> 404 Four Sensor Alert &amp; 405 Two Sensor Alert, collision warning system, (auditory)</p> <p><b>LUCI:</b> Incline/tip warning (auditory)</p>
					Backup camera gives visual display of potential collisions or other hazards prior to driving in reverse	<p><b>Cheelcare:</b> Aware A1 / A2 / A3 backup Cameras</p> <p><b>Tadibrothers:</b> backup camera</p> <p><b>Quantum Rehab:</b> backup camera</p>

✓ Human driver monitors the driving environment

# Let's Talk Sensors

- Variety of sensors...but each type has its strengths and weaknesses!
- Examples:
  - Stereo Vision Cameras
  - Radar
  - Ultrasonic
  - LIDAR



# Warning Systems: No Automation



## Adaptive Switch Labs

- ✓ Photo-electric switches in a 2 or 4 switch array, provides auditory feedback when approaching an obstacle
- ✓ Can be programmed to 'alert' between 4-21" from obstacle (user reaction time)
- ✓ Driver must hear and modify driving based on warning
- ✓ Does **not** detect drop-offs



Close Up View

# Warning Systems: No Automation



- ✓ Ultrasonic blind spot sensors: mounted at user's preferred location on PWC
- ✓ Feedback choice of visual, auditory, and/or vibration (up to 3 pads)
- ✓ Most robust model, Sentina, provides 180 degrees horizontal/rearview and 50 degrees vertical coverage
- ✓ Can add additional blind spot coverage via Echo Head (up to 3)
- ✓ Does **not** detect drop-offs or soft material

SENTINA ULTRASONIC SENSOR UNIT



ECHO HEAD



VIBRATORY PAD

VISUAL/  
AUDITORY  
FEEDBACK  
DISPLAY PANEL





# Braze Mobility



# Warning Systems: No Automation

**LUCI** : Incline/tip warning

- ✓ Audible alert when LUCI detects the chair is driving at an unsafe angle where the chair is at risk for tipping over



# Warning Systems: No Automation

## Backup Cameras



“Last week, I received the Quantum<sup>®</sup> wheelchair backup camera. I honestly didn't know what to expect, but in the first few days of having it, I'm pleasantly surprised. I don't enjoy backing up my Quantum Rehab Wheelchair, but I'm finding now that I do it a lot because I officially have eyes in the back of my head! What's your superpower?”

“Having the power wheelchair backup camera has given me peace of mind, saved my walls and other people's toes! I have a service dog that I use both at home and in public. With the backup camera, I can safely find his paws and not accidentally hurt him. I'm sure he appreciates that!”

# SMART WHEELCHAIR (SWC) CONTINUUM

Human driver and automated driving system monitor the driving environment						
1	<b>Driver Assistance:</b> Single function assisted navigation (Speed or Steering)	An active system that can make adjustments to inputs for only one function (e.g., either speed or steering) to assist with navigation. The system DOES intervene.	Collaborator-the user is in control with assistance from the system	Driver controls all driving options except for emergency stops in response to detected collisions or other hazards	Sensors that warn driver of potential collisions or other hazards and the system stops the wheelchair if the driver does not respond appropriately	NA
				Increases driving efficiency by reducing compensatory movement, for example driving on a side slope.	Tracking technology: reduces joystick movements or switch activations and reduces time to move between locations	<b>Quantum Rehab:</b> <u>Accu-Trac</u> <b>Invacare:</b> G-Trac <b>Permobil:</b> <u>ESP</u> <b>Sunrise Medical:</b> <u>SureTrac</u> <b>AMYLOR:</b> Smart-Track

- ✓ Level 1: Human driver and automated driving system monitors the driving environment

# Driver Assistance: Emergency Stop

permobil

Remote Stop



Remote Emergency Stop Switch



# Driver Assistance: Tracking Technology

- ✓ Increased efficiency! Reduces joystick movement/switch activations required to get from Point A to Point B
- ✓ Optional (meaning it can be denied) and is often not ordered/approved for clients who can benefit

**permobil**

Enhanced Steering Performance (ESP)

**INVACARE**

G-Trac

**SUNRISE MEDICAL**

SureTrac

**AMYLOR**

Smart-Track

**QUANTUM REHAB**

Accu-Trac



# SMART WHEELCHAIR (SWC) CONTINUUM

SWC Level	Name	Definition	Human Role	Functional Example	Product Feature Example	Product Example
<b>Human driver and automated driving system monitor the driving environment, cont.</b>						
2	<b>Advanced Driver Assistance:</b> Multiple function assisted navigation. (Speed and Steering)	An active system that can make both speed and steering adjustments simultaneously to the driver's inputs to assist with navigation. The system DOES intervene.	Cooperator- The user monitors and engages while the system can adjust inputs.	Driver can steer. System will avoid collisions, drop-offs, and/or tipping by Simultaneously controlling speed and direction	Driver can continue driving, but not in the direction of a hazard. System imposes a restriction in travel that can be overridden	<b>LUCI:</b> Navigation assistance/collision avoidance, Drop-off protection
				Driver can control speed. System will automatically slow, as needed	Driver can increase speed, but system will slow in response to environment, such as walking in a crowd	<b>LUCI:</b> Crowd confidence and dynamic slowing

- ✓ Level 1: Human driver and automated driving system monitors the driving environment

# Advanced Driver Assistance

**LUCI** : Fusion Sensors offer a 360-degree view of the driving environment detecting obstacles AND drop-offs


- **Light Blue** represents coverage by stereo vision cameras and infrared
- **Purple** represents radar coverage
- **Green** represents ultrasonic coverage



# Advanced Driver Assistance

## LUCI : Steering + Speed:

- Sensors will slow the chair as it approaches an obstacle giving the driver an opportunity to self correct/make a course deviation
- If the driver does not respond (keeps the same input) or is not able to respond in time, LUCI will intervene and stop the chair






# Advanced Driver Assistance

## LUCI : Drop-off Detection


- If LUCI's FusionSensors detect an unsafe **drop off**, speed will automatically slow and then stop, if needed, to prevent dropping off a curb or the edge of a ramp



# REACTION TIME FEATURE

Begins Deceleration Sooner

SLOWER REACTION TIME (>1.0 SECONDS)



SAME STOPPING POINT



Begins Deceleration Later

FASTER REACTION TIME (>0.5 SECONDS)



# Advanced Driver Assistance

**LUCI** : Crowd Confidence



# Advanced Driver Assistance

LUCI :White Cane Filter



## Advanced Driver Assistance

LUCI :White Cane Filter



# SMART WHEELCHAIR (SWC) CONTINUUM

Automated driving system monitors the driving environment						
3	<b>Conditional Automation:</b> Autonomously navigate through a specific process and adapt under specific conditions.	An active system that makes limited, fully automated actions in response to the user inputs. The system DOES intervene.	Initiator/Supervisor- Users must be ready to drive when autonomous features are not engaged.	Ability to navigate to a destination. The driver can initiate and stop movement, as desired, but stopping is not required	System follows a preprogrammed 'map' or tape on the floor and modifies driving in response to sensor feedback	<b>Smile Smart System (SSS):</b> driver initiates and stops movement with <u>switch</u>  <b>LUCI:</b> <u>RampAssist™</u>

# Conditional Automation

## SMILE Smart System

- ✓ Driver initiates and stops movement with switch
- ✓ The PWC follows a tape track that can be used indoors and/or outdoors
- ✓ Sensors prevent collision (anti-collision sensors)
- ✓ Line following can be used as a safe pathway from which to develop switch access skills and joystick use over time, gently adapting settings as personal abilities evolve



# Conditional Automation



# SMART WHEELCHAIR (SWC) CONTINUUM

SWC Level	Name	Definition	Human Role	Functional Example	Product Feature Example	Product Example
4	<b>Highly Autonomous System:</b> The wheelchair is fully autonomous for an entire trip in specific environments.	An active system where driver input is unnecessary in specific environments and situations.	Occupant in specific environments-no human interaction needed	Ability to navigate to a destination while deciding an optimal process for negotiating obstacles and terrain. The system controls all features in specific environments	The system controls all features in specific environments	<b>NA</b>
5	<b>Fully Autonomous System:</b> The wheelchair can navigate without a human in all environments.	An active system where driver input is not required.	Occupant in all environments -no human interaction needed	Autonomous in all environments -The system controls all features, everywhere, <u>at all times</u> , in all conditions.	The system controls all features	<b>NA</b>

# Highly/Fully Autonomous System

No FDA approved medical devices currently on the market!

But...

# Connectivity

A SMART wheelchair also means the ability to integrate to the connected world


# Connectivity: BT to other devices

All Complex Rehab power wheelchairs offer Bluetooth (BT) that allows switch access or mouse emulation to external devices, i.e computer, tablet, smartphone, SGD

- ✓ Invacare: LiNX
- ✓ Permobil: R-net
- ✓ Quickie: R-net
- ✓ Quantum: Q-Logic
- ✓ AmyLior: R-net
- ✓ Merits: R-net




# Connectivity: BT to Monitoring Apps



permobil

## MyPermobil


- Real-time battery status, distance traveled, seating activity
- Integrated map with GPS location
- Voice Assistant (e.g. Alexa)
- Fleet Management



SUNRISE  
MEDICAL.

## Switch-It

- Monitors time spent in various seating positions + alerts when time for position change
- Share with care team to create individualized pressure relief programs




LUCI

## MyLUCI App

- View collisions, drop-offs avoided
- View live sensor feedback with LUCIView™
- Manage health and seating alerts
- Invite caregivers, share data, trigger event notifications
- Voice Assistant (e.g. Alexa)
- Driver location (using GPS)
- Battery usage

Finally, a chair with a LUCI View™



# Connectivity: Health and Wellness

## MyLUCI App:


- Apple Health and Google Fit compatible heart rate monitoring
- Heart rate tracking and notification of elevated heart rate sent to care team



# Connectivity: WiFi

## LUCI:

- ✓ Over-the-air updates are pushed automatically- LUCI updates overnight to the latest software version
  - New product features added via software update, i.e. RampAssist™
- ✓ Tech support can “see” what LUCI “sees” to assist with troubleshooting.



## Permobil:

- ✓ QuickConfig
  - Allows programming and customization of the chair, i.e. memory seating positions, standing sequence, drive profiles




# Judging Smart White Paper

- White paper designed around questions that users should be asking of their assistive technology
- This has also been adapted into a downloadable slide deck, perfect for clinicians/ATPs/users wanting to boost their language and understanding of SMART Assistive Tech.

Download both at:

[www.luci.com/smart/](http://www.luci.com/smart/)



# The Person

Clinical Applications

# Clinical Applications: Accessibility

- Maneuvering around tight spaces is very difficult
  - Examples:
    - up a ramp to get into a van
    - within a van to line up with tie downs or a locking system
    - a crowded hallway at school or an assisted living facility
    - airport
    - grocery store



# Clinical Applications: Accessibility



# Clinical Applications: Obstacles

- To avoid obstacles, the driver must see them
  - Hard to see areas: behind and low
- Also need to gauge distance and respond in a timely manner
- These are often very difficult for our PWC users



# Clinical Applications: Motor, Visual, and Cognitive Requirements

- Motor limitations may limit driving precision and reaction time
- Visual limitations may make driving more difficult, specifically lack of acuity and visual spatial concerns (i.e. depth perception), visual field cut
- Visual field neglect or inattention
- Cognitive limitations may lead to a reduced understanding of the implications of certain driving maneuvers, such as driving off of a curb or colliding with an obstacle
- **Many clients have more than one area of involvement**



“The envelope of who is going to be able to safely operate a wheelchair in a whole variety of environments just opened right up”

JEAN MINKEL, PT/ATP  
SENIOR VICE PRESIDENT  
INDEPENDENCE CARE SYSTEMS  
NEW YORK, NY

# Clinical Applications: Summary

- So, who can benefit from Smart Wheelchair technologies?
  - Anyone who is not driving efficiently and safely, to their full potential
  - A client who has been deemed unsafe to drive a PWC and is currently in a dependent situation, i.e. tilt-in-space
  - A client who requires assistance/intermittent assistance with management of a MWC or a PWC

NOTE: Important to match client needs with specific product parameters

- There is no one product that meets everyone's needs, that is why understanding the continuum is important

---

Let's Try IT!

# Questions?

## Other Resources:

- 1.Scott Crawford: [Orientation & Mobility for Visually Impaired Wheelchair Users](#)
- 2.American Printing House (APH) for the Blind: [APH: Orientation & Mobility for Wheelchair Users](#)
- 3.Textbook: [Foundations of Orientation and Mobility, Volume I: History and Theory | American Printing House](#)

# LUCI Contact Information

- Lindsey Sharpe, PT, DPT, ATP: Regional Director of Clinical Education
  - Email: [lindsey@luci.com](mailto:lindsey@luci.com)
- Manfred Call VI: Territory Manager
  - Email: [call6manfred@gmail.com](mailto:call6manfred@gmail.com)